

Specialty Syringes and Accessories



Threaded Plunger Syringes

For applications requiring extremely precise plunger movement LT Terminus Volumes from 25 µl to 10 ml

Threaded Plunger Syringes			
Catalog No.	Volume	Terminus	Dispense Volume
CGS 8446.63	25 µl	LT	0.33 µl/revolution
CGS 8447.63	50 µl	LT	0.66 µl/revolution
CGS 8448.63	100 µl	LT	1.32 µl/revolution
CGS 8449.63	250 µl	LT	3.31 µl/revolution
CGS 8450.63	500 µl	LT	6.62 µl/revolution
CGS 8451.63	500 µl	LT	5.29 µl/revolution
CGS 8452.63	1 ml	LT	13.23 µl/revolution
CGS 8453.63	2.5 ml	LT	37.79 µl/revolution
CGS 8454.63	5 ml	LT	75.57 µl/revolution
CGS 8455.63	10 ml	LT	151.14 µl/revolution

Popper Multifit Glass Syringes

These are standard glass Syringes with Luer Lock tips. The barrels and plungers of the same-sized Syringes are interchangeable.

Popper Multifit Glass Syringes			
Catalog No.	Syringe Size	Catalog No.	Syringe Size
CGS 8456.63	2 cc	CGS8459.20	30 cc
CGS 8457.63	20 cc	CGS8460.20	50 cc
CGS 8458.63	5 cc	CGS 8461.20	10 cc

Becton-Dickinson Yale Glass Syringe with Robb Tip, 100 ml

The plunger and barrel of this glass Syringe are ground and mated for a tight fit. They are not interchangeable with other barrels or plungers. The parts are numbered for easy reassembly.

This Syringe dispenses large volumes in short time periods. Therefore, it has a Robb tip with a 12 gauge bore of 0.223 cm (0.088 in). Fitted with Luer lock tip. 100 ml size only.

Catalog No.	Product
CGS 8479.63	Becton-Dickinson Yale Glass Syringe with Robb Tip, 100 ml



Syringe Guides				
Catalog No.	Volume	For Use with Syringe Series		
		700	1700	7000
CGS 8462.63	0.5 to 5 µl			•
CGS 8463.63	5 to 10 µl	•		
CGS 8464.63	5 to 10 µl		•	
CGS 8465.63	25 to 500 µl	•		
CGS 8466.63	25 to 500 µl		•	



Needle Cleaning Kit

Fine gauge tungsten wires used to clear plugged needles. Biodegradable cleaning solutions are used to remove residues from syringe needles and barrels. All cleaning wires are 7 inches long and are

sold in a package of 10. Cleaning concentrate is sold separately.

Cleaning Wires		
Catalog No.	To Clean Needle Gauge	Wire Size OD
CGS 8467.63	23s, 26s, 31 to 33	0.0762 mm (0.00300 in)
CGS 8468.63	26s, 31 to 33	0.0889 mm (0.00350 in)
CGS 8469.63	22s, 25s, 28 to 30	0.1262 mm (0.00497 in)
CGS 8470.63	27	0.1674 mm (0.00659 in)
CGS 8471.63	24 to 26	0.2070 mm (0.00815 in)
CGS 8472.63	22, 23	0.3066 mm (0.01207 in)

Catalog No.	Product
CGS 8473.63	Cleaning Concentrate, 70 ml
CGS 8474.63	Cleaning Concentrate, 500 ml
CGS 8475.63	All Cleaning Wires and 70 ml of Cleaning Concentrate



Removable (RN) Adapters

Provides connection of Luer fittings to RN syringes ≥ 250 µl

Catalog No.	Product
CGS 8476.63	RN Hub to Female Luer
CGS 8477.63	RN Hub to Male Luer
CGS 8478.63	RN Hub to Male Luer Lock