

## EasyMount Vertical Diffusion Chamber System

EasyMount Diffusion Chambers and Systems			Inserts for EasyMount Diffusion Chambers				
Catalog No.	Model	Product	Catalog No.	Model	Aperture	Area	Description
<b>Diffusion Chambers</b>			<b>For P2300 EasyMount Diffusion Chambers</b>				
CGS 8923.75	P2300	EasyMount Diffusion Chamber	CGS 8940.75	P2302	12.7 mm dia	1.13 cm <sup>2</sup>	For Snapwell cell culture
CGS 8924.75	P2250	EasyMount Diffusion Chamber, 2 to 8 ml	CGS 8941.75	P2304	2.8 x 11 mm	0.3 cm <sup>2</sup>	For mouse intestine
CGS 8925.75	P2400	EasyMount Low Volume Diffusion Chamber, ~1 ml	CGS 8942.75	P2305	4.5 x 11.2 mm	0.5 cm <sup>2</sup>	For rat and larger animal intestine
CGS 8926.75	P2500	EasyMount Diffusion Chamber for Cell Culture	CGS 8943.75	P2306	2.8 x 1.5 mm	0.04 cm <sup>2</sup>	For mouse trachea, similar tissues
<b>Diffusion Chamber Systems</b>			CGS 8944.75	P2307	2 mm dia	0.031 cm <sup>2</sup>	For small tissues (i.e., biopsies)
CGS 8927.75	EM-CSYS	EasyMount Diffusion Chamber System with Snapwell Cell Culture Inserts	CGS 8945.75	P2310	5 mm dia	0.2 cm <sup>2</sup>	With mounting pins
CGS 8928.75	EM-RSYS	EasyMount Diffusion Chamber System for Intestine	CGS 8946.75	P2311	6.2 mm dia	0.3 cm <sup>2</sup>	With mounting pins
CGS 8929.75	EM-BSYS	EasyMount Biopsy, Small Tissue Diffusion Chamber System	CGS 8947.75	P2312	8 mm dia	0.5 cm <sup>2</sup>	With mounting pins
CGS 8930.75	EM-LVSY	EasyMount Low Volume Diffusion Chamber System	CGS 8948.75	P2313	9.5 mm dia	0.71 cm <sup>2</sup>	With mounting pins
CGS 8931.75	EM-CELL	EasyMount Diffusion Chamber System for Cell Culture	CGS 8949.75	P2314	11.3 mm dia	1 cm <sup>2</sup>	With mounting pins
<b>Accessories for Diffusion Systems</b>			CGS 8950.75	P2315	12.7 mm dia	1.26 cm <sup>2</sup>	With mounting pins
CGS 8932.75	P2000	EasyMount Stand including base, legs, heat block, cradles, cam levers, air valves, tubing and fittings	<b>For P2250 EasyMount Diffusion Chamber, 2 to 8 ml</b>				
CGS 8933.75	P2020	Voltage and Current Electrode Set	CGS 8951.75	P2252	5 x 22 mm	1 cm <sup>2</sup>	For tissues, ie: rat, guinea pig or rabbit intestine
CGS 8934.75	P2020-V	Voltage Electrode Set	<b>For P2400 EasyMount Diffusion Chambers, Low Volume</b>				
CGS 8935.75	P2020-I	Current Electrode Set	CGS 8952.75	P2404	2.8 x 9.2 mm	0.25 cm <sup>2</sup>	For mouse intestine
CGS 8936.75	P2023-20	Electrode Tips, pkg. of 20	CGS 8953.75	P2405	4.5 x 9.2 mm	0.4 cm <sup>2</sup>	For rat and larger animal intestine
CGS 8937.75	P2023-50	Electrode Tips, pkg. of 50	CGS 8954.75	P2406	2.8 x 1.5 mm	0.04 cm <sup>2</sup>	For mouse trachea, similar tissues
CGS 8938.75	P2023-100	Electrode Tips, pkg. of 100	CGS 8955.75	P2407	2 mm dia	0.031 cm <sup>2</sup>	For small tissues (i.e., biopsies)
CGS 8939.75	P2024	Electrodes Leads, 17.8 cm (7 in), pkg. of 4	CGS 8956.75	P2410	5 mm dia	0.2 cm <sup>2</sup>	With mounting pins
			CGS 8957.75	P2411	6.2 mm dia	0.3 cm <sup>2</sup>	With mounting pins
			CGS 8958.75	P2412	8 mm dia	0.5 cm <sup>2</sup>	With mounting pins
			CGS 8959.75	P2413	9.5 mm dia	0.71 cm <sup>2</sup>	With mounting pins
			<b>For P2500 EasyMount Diffusion Chambers for Cell Culture</b>				
			CGS 8960.75	P2504	2.8 x 9.2 mm	0.25 cm <sup>2</sup>	For mouse intestine
			CGS 8961.75	P2505	4.5 x 9.2 mm	0.4 cm <sup>2</sup>	For rat and larger animal intestine
			CGS 8962.75	P2524	6.5 mm	0.33 cm <sup>2</sup>	For Transwell™ cell culture
			CGS 8963.75	P2525	12 mm	1.13 cm <sup>2</sup>	For Transwell™ cell culture
			CGS 8964.75	P2530	6 mm	0.3 cm <sup>2</sup>	For Falcon cell culture
			CGS 8965.75	P2531	10 mm	0.9 cm <sup>2</sup>	For Falcon cell culture
			CGS 8966.75	P2535	8 mm OD	0.5 cm <sup>2</sup>	For Nunc Anopore cell culture
			CGS 8967.75	P2538	12 mm OD	0.6 cm <sup>2</sup>	For Millicell cell culture