


Sealpette Electronic Single Channel Pipettors



Sealpette Electronic Single Channel Pipettor

- Microprocessor control eliminates operator variability
- Reduces risk of RSI (repetitive strain injury)
- Clear LCD display shows current function, volume and piston speed
- One lightweight instrument covers all the essential liquid handling operations: pipetting, reverse pipetting, multiple dispensing (for sample replicates and reagent dispensing), dilute dispensing (up to 1:24 dilution) and mixing
- Tip ejector precisely placed for single-handed use

These lightweight electronic Sealpettes look and feel like manual pipettors, yet they provide the advanced liquid handling features that only microprocessor control can provide. All features are accessed through a simple six button keypad.

Charging Stands

A charging stand is an essential part of the electronic Sealpette pipettor. Two stands are available. An individual charging stand to hold one Sealpette, or a 4-position carousel stand. The Sealpette electronic pipettor is available with or without the single charging stand. Should you prefer the 4-position carousel stand, simply purchase a Sealpette without a stand and order one of the 4-Position carousel charging stands to go with it.

Electronic Sealpettes Single Channel

110 VAC, 60 Hz	240 VAC, 50 Hz, UK Plug	Other Power*	Volume Range	Increment	Inaccuracy	Imprecision
With Charging Stand						
CGS 8866.33	CGS-8878.33	CGS-8900.33	0.2 to 10 µl	.1 µl	0.90%	0.50%
CGS 8867.33	CGS-8879.33	CGS-8901.33	5 to 100 µl	µl 1	0.40%	0.15%
CGS 8868.33	CGS-8880.33	CGS-8902.33	10 to 250 µl	µl 1	0.40%	0.15%
CGS 8869.33	CGS-8881.33	CGS-8903.33	20 to 500 µl	µl 5	0.40%	0.15%
CGS 8870.33	CGS-8882.33	CGS-8904.33	50 to 1000 µl	µl 5	0.40%	0.15%
CGS 8871.33	CGS-8883.33	CGS-8905.33	100 to 5000 µl	10 µl	0.50%	0.15%
Without Charging Stand						
CGS 8872.33	CGS-8884.33	CGS-8906.33	0.2 to 10 µl	.1 µl	0.90%	0.50%
CGS 8873.33	CGS-8885.33	CGS-8907.33	5 to 100 µl	µl 1	0.40%	0.15%
CGS 8874.33	CGS-8886.33	CGS-8908.33	10 to 250 µl	µl 1	0.40%	0.15%
CGS 8875.33	CGS-8887.33	CGS-8909.33	20 to 500 µl	µl 5	0.40%	0.15%
CGS 8876.33	CGS-8888.33	CGS-8910.33	50 to 1000 µl	µl 5	0.40%	0.15%
CGS 8877.33	CGS-8889.33	CGS-8911.33	100 to 5000 µl	10 µl	0.50%	0.15%

* Note: Must specify desired power/plug configuration.



Single Pipettor Charging Stand



4-Position Carousel Charging Stand

Catalog No.

Product

- CGS 8912.33** Single Pipettor Charging Stand, 110 VAC, 60 Hz
- CGS 8913.33** Single Pipettor Charging Stand, 240 VAC, 50 Hz UK Plug
- CGS 8914.33** Single Pipettor Charging Stand, other power, must specify desired power/plug configuration
- CGS 8915.33** 4-Position Carousel Charging Stand, 110 VAC, 60 Hz
- CGS 8916.33** 4-Position Carousel Charging Stand, 240 VAC, 50 Hz UK Plug
- CGS 8917.33** 4-Position Carousel Charging Stand, other power, must specify desired power/plug configuration