


## Sealpette Tip Selection Guide



Sealtip Selector for Manual Sealpettes			
Catalog No.	Pipettor	Range	Suitable Tips
CGS 8981.33	1 Channel Variable Micro Vol.	0.5 to 10 µl	A, D, E, J, H
CGS 8982.33	1 Channel Variable Volume	2 to 20 µl	A, E, J
CGS 8983.33	1 Channel Variable Volume	5 to 50 µl	A, E, J
CGS 8984.33	1 Channel Variable Volume	10 to 100 µl	A, E, J
CGS 8985.33	1 Channel Variable Volume	20 to 200 µl	A, E, J
CGS 8986.33	1 Channel Variable Volume	50 to 200 µl	A, E, J
CGS 8987.33	1 Channel Variable Volume	100 to 1000 µl	BF, K,
CGS 8988.33	1 Channel Variable Volume	200 to 1000 µl	BF, K,
CGS 8989.33	1 Channel Variable Volume	1000 to 5000 µl	CG,
CGS 8990.33	4 Channel Variable Volume	5 to 50 µl	EJ, M
CGS 8991.33	4 Channel Variable Volume	50 to 250 µl	EJ, M
CGS 8992.33	8 Channel Variable Micro Vol.	0.5 to 10 µl	DH
CGS 8993.33	8 Channel Variable Volume	5 to 50 µl	EJ, M
CGS 8994.33	8 Channel Variable Volume	50 to 250 µl	EJ, M
CGS 8995.33	12 Channel Variable Volume	0.5 to 10 µl	DH
CGS 8996.33	12 Channel Variable Volume	5 to 50 µl	EJ, M
CGS 8997.33	12 Channel Variable Volume	50 to 250 µl	EJ, M
CGS 8998.33	1 Channel Fixed Volume	5 µl	A, D, E, J, H
CGS 8999.33	1 Channel Fixed Volume	10 µl	A, D, E, J, H
CGS 8000.39	1 Channel Fixed Volume	20 µl	A, E, J
CGS 8001.39	1 Channel Fixed Volume	25 µl	A, E, J
CGS 8002.39	1 Channel Fixed Volume	50 µl	A, E, J
CGS 8003.39	1 Channel Fixed Volume	100 µl	A, E, J
CGS 8004.39	1 Channel Fixed Volume	200 µl	A, E, J
CGS 8005.39	1 Channel Fixed Volume	250 µl	A, E, J
CGS 8006.39	1 Channel Fixed Volume	500 µl	BF, K
CGS 8007.39	1 Channel Fixed Volume	1000 µl	BF, K
CGS 8008.39	1 Channel Fixed Volume	2000 µl	CG
CGS 8009.39	1 Channel Fixed Volume	5000 µl	CG

Sealpette Tips					
Catalog No.	Tip Ref	Color	Pack Size	Name	Range
CGS 8968.33	A	Yellow	1000 pack	Economy	0.5 to 250 µl
CGS 8969.33	B	Blue	100 pack	Economy	200 to 1000 µl
CGS 8970.33	C	White	300 pack	Economy	1000 to 5000 µl
CGS 8971.33	D	White	1000 pack	Ultra-micro	0.5 to 10 µl
CGS 8972.33	E	White	1000 pack	Superfine	0.5 to 250 µl
CGS 8973.33	F	White	400 pack	Superfine	200 to 1000 µl
CGS 8974.33	G	White	100 pack	Superfine	1000 to 5000 µl
CGS 8975.33	H	White	10x96 rack	Ultra-micro	0.5 to 10 µl
CGS 8976.33	J	White	10x96 rack	Superfine	0.5 to 250 µl
CGS 8977.33	K	White	10x96 rack	Superfine	200 to 1000 µl
CGS 8978.33	L	White	10x96 rack	Superfine	50 to 1200 µl
CGS 8979.33	L	White	1000 pack	Superfine	50 to 1200 µl
CGS 8980.33	M	White	100 pack	4 band	0.5 to 250 µl

Sealtip Selector for Electronic Sealpettes			
Catalog No.	Pipettor	Range	Suitable Tips
CGS 8010.39	1 Channel Electronic	0.2 to 10 µl	D H
CGS 8011.39	1 Channel Electronic	5 to 100 µl	E J
CGS 8012.39	1 Channel Electronic	10 to 250 µl	E J
CGS 8013.39	1 Channel Electronic	20 to 500 µl	F K
CGS 8014.39	1 Channel Electronic	50 to 1000 µl	F K
CGS 8015.39	1 Channel Electronic	100 to 500 µl	C,
CGS 8016.39	4 Channel Electronic	5 to 100 µl	E, J, M
CGS 8017.39	4 Channel Electronic	25 to 250 µl	E, J, M
CGS 8018.39	8 Channel Electronic	0.2 to 10 µl	D H
CGS 8019.39	8 Channel Electronic	5 to 100 µl	E, J, M
CGS 8020.39	8 Channel Electronic	25 to 250 µl	E, J, M
CGS 8021.39	8 Channel Electronic	50 to 1200 µl	L
CGS 8022.39	12 Channel Electronic	0.2 to 10 µl	D H
CGS 8023.39	12 Channel Electronic	5 to 100 µl	E, J, M
CGS 8024.39	12 Channel Electronic	25 to 250 µl	E, J, M