



Electrophoresis Systems

CoZap™



Using CoZap is easy:

1. Place the CoZap pad in your destaining tank.
2. Remove the gel after destaining.



BEFORE

AFTER

CoZap

CoZap is used for rapid removal of Coomassie Blue stain from electrophoresis gels without the need to change the destaining solution. CoZap is a unique pad that has a high absorbance for Coomassie blue stain and is thus very effective in destaining gels. CoZap absorbs any free dye in the solution making gel destaining 20% faster than with conventional method. It is one of the most effective destaining methods on the market.

- No need to change the destaining solution!
- Fast and simple
- No charcoal or dye residues
- No subsequent destaining required
- One pad can destain up to 10 gels
- 20% faster than conventional methods

Catalog No.	Product	Quantity
CGS 8395.39	CoZap Pads	25
CGS 8396.39	CoZap Pads	100