


Ultrasonic Baths



- Baskets prevent objects touching the sides of the tank and allow items to be placed in or removed from the bath without immersing your hands in the chemical solution
- Lids minimize operating noise and reduce evaporation
- The use of a basket is recommended

Suitable Lids and Baskets for UB Ultrasonic Baths

Bath Model	Suitable Lid	Suitable Basket	Basket Dimensions, L x W x D	Overall Depth to Tank Rim
2 Liter	CGS-8418.78	CGS-8423.78	10.2 x 8.9 x 6.4 cm (4 x 3.5 x 2.5 in)	9.3 cm (3.7 in)
3 Liter	CGS-8419.78	CGS-8424.78	19.5 x 12 x 6.4 cm (7.6 x 4.7 x 2.5 in)	9.3 cm (3.7 in)
6 Liter	CGS-8420.78	CGS-8425.78	23.7 x 10.2 x 10.2 cm (9.3 x 4.0 x 4.0 in)	13 cm (5.1 in)
14 Liter	CGS-8421.78	CGS-8426.78	26.3 x 23.6 x 10.2 cm (10.4 x 9.3 x 4.0 in)	13 cm (5.1 in)
22 Liter	CGS-8422.78	CGS-8427.78	44.5 x 25.2 x 10.2 cm (17.5 x 9.9 x 4.0 in)	13 cm (5.1 in)

Ultrasonic Baths

	2 Liter	3 Liter	6 Liter	14 Liter	22 Liter
Overall Dimensions, L x W x D	17.6 x 16.2 x 25 cm (6.9 x 6.4 x 9.8 in)	26.5 x 16.2 x 25 cm (10.4 x 6.4 x 9.8 in)	20.5 x 32.5 x 31 cm (8.1 x 12.8 x 12.2 in)	38 x 32.5 x 31 cm (15 x 12.8 x 12.2 in)	55.5 x 32.5 x 31 cm (21.9 x 12.8 x 12.2 in)
Working Volume, L x W x D	14.9 x 13.5 x 10 cm (5.9 x 5.3 x 3.9 in)	23.8 x 13.5 x 10 cm (9.4 x 5.3 x 3.9 in)	15 x 30 x 14 cm (5.9 x 11.8 x 5.5 in)	32.5 x 30 x 14 cm (12.8 x 11.8 x 5.5 in)	50.5 x 30 x 14 cm (19.9 x 11.8 x 5.5 in)
Drain Facility	—	—	—	Drain tap	Drain tap
Typical Temperature, Switchable	—	—	50°C	50°C	50°C
Heater Power	—	—	200 W	600 W	800 W
Transducers	1	2	3	6	9
Generators	1	1	1	2	3
Power RMS:	63 W	75 W	89 W	162 W	238 W
Peak	125 W	150 W	180 W	325 W	475 W
Overall Consumption	60 W	520 W	300 W	800 W	1100 W
Typical Noise Level A Leq* dB	< 65	< 68	< 68	< 70	< 72
Variable Timer (Sonics)	—	—	0 to 30 mins	0 to 30 mins	0 to 30 mins
Weight	3.8 kg (8.4 lb)	5.8 kg (12.8lb)	6.3 kg (13.9 lb)	11.1 kg (24.4 lb)	15.7 kg (34.5 lb)
Safety Temperature	—	—	Fixed cut-out	Fixed cut-out	Fixed cut-out
Catalog No.	CGS-8428.78	CGS-8429.78	CGS-8430.78	CGS-8431.78	CGS-8432.78