



## **CGS 4004.20**

### **Manual or all-automatic solid sample feed**

Whether manual or fully automated, the solid samplers allow direct introduction of solid samples into the graphite furnace. Uniquely, fast and efficiently – this system opens up an entirely new field of AAS analyses.

#### Manual sample feed – SSA6

- Manually guided, reproducible introduction of the platform
- Automatic transfer of weight data to software

#### Automatic sample feed – SSA61

- Fully automatic introduction of loaded platforms into the graphite furnace
- Integrated microbalance for automated weight data transfer
- Different operating modes for up to 10 sample places available
- Required for the use of the solids technique in routine mode