

Mistura de Solvente

Temperature: 40° (2.0 min. hold) (6°/min.) -225°C
Injector: 225°C
Detector: MSD
Carrier Gas: 28 cm/sec., He

- | | |
|-------------------------|----------------------------|
| 1- pentane | 14- trichloroethylene |
| 2- hexane | 15- acetonitrile |
| 3- MTBE | 16- tetrachloroethylene |
| 4- heptane and methanol | 17- chloroform |
| 5- cyclohexane | 18- 2-butanol |
| 6- acetone | 19- toluene |
| 7- THF | 20- propanol |
| 8- trichloroethane | 21- 1,4-dioxane |
| 9- ethyl acetate | 22- isobutyl alcohol |
| 10- MEK and ethanol | 23- pyridine |
| 11- methylene chloride | 24- n,n-dimethyl formamide |
| 12- benzene | 25- bromoform |
| 13- butyl cellosolve | |

CGS 1001.18

