

Operating Tables and Accessories

Plastic Operating Table



This table is molded entirely of plastic and is for use with animals weighing up to 4.5 kg (10 lb).

The 24.5 x 37 cm (9-3/4 x 14-3/4) table has a 6.4 mm (1/4 in) raised central platform, 6.4 x 20.3 cm (2-1/2 x 8 in), to support the animal's body and head.

Sloping drainage channels on each side collect body fluids and deliver them to the end of the table.

The table is 5.3 cm (2-1/8 in) thick. Folding legs increase the height to 9.5 cm (3-3/4 in). The legs fold independently to provide versatility in positioning the animal at an optimum angle.

The table is supplied with 5.5 meters (18 ft) of elastic cord and 20 rings. The cord and rings are made of smooth, white, ABS plastic that is easy-to-clean, corrosion-proof and cannot be penetrated by radioactive substances. The corners and sides of the table are slotted to receive elastic cord and ring assemblies used to tie down the animal.

Catalog No.	Product
CGS 8006.23	Plastic Operating Table

Frog Board




To prevent warping, this 9 x 14 cm (3-1/2 x 5-1/2 in) Frog Board is made of a thick piece of cork cemented to a heavy piece of 12.5 mm (1/2 in) thick white acrylic. A clamp and locking screw are mounted on the side with an opening to fit a standard 9.5 mm (3/8 in) diameter rod

(not supplied).

Catalog No.	Product
CGS 8007.66	Frog Board

Stainless Steel Thoracic Positioner



This Thoracic Positioner will hold the animal in the desired operating position. Made of heavy gauge stainless steel for durability, this positioner has sure grip rubber feet to hold it in place on any flat-top surface or operating table. It measures, H x W x D, 8.25 x 30.5 x 26.7 cm (3-1/4 x 12 x 10-1/2 in) and weighs 2.25 kg (5 lb).

Catalog No	Product
CGS 8008.66	Stainless Steel Thoracic Positioner

Plastic Thoracic Positioners



- Durable, lightweight plastic
- Transparent to light and x-rays
- Easy-to-clean
- Low heat loss

These Thoracic Positioners are designed to hold the animal in the desired operating position and are available in three sizes.

Since these positioners are constructed of plastic instead of stainless steel, the animal will lose less body heat during surgery.

Plastic Thoracic Positioners			
Catalog No.	For use with	Animal Weight	Dimensions, H x W x D
CGS 8009.66	Cats and Small Dogs	Up to 9 kg (20 lb)	5.7 x 20.3 x 17.8 cm (2-1/4 x 8 x 7 in)
CGS 8010.66	Medium Dogs	9 to 18 kg (20 to 40 lb)	7 x 25.4 x 20.5 cm (2-3/4 x 10 x 8 in)
CGS 8011.66	Large Dogs	18 to 27.2 kg (40 to 60 lb)	7.6 x 35.6 x 22.7 cm (3 x 14 x 9 in)

CareFRESH® Absorbent Sheets



- Super absorbent
- Inexpensive
- Ideal for post surgery, maintaining body heat during recovery, protecting exam tables, and for cleanup of accidental spills

These disposable sheets are ideal for numerous applications. They help prevent body heat loss after surgery. When used under small animal cages, they help make excrement cleanup fast and easy. They can also be used to protect examination tables. Made from a super absorbent cellulose fiber, they are also ideal for cleaning up various spills and animal accidents.

An individual sheet measures 40.6 x 50.8 cm (16 x 20 in). Select from 6.4 mm (1/4 in) or 3.2 mm (1/8 in) thickness. They are supplied in a plastic bag.

Catalog No.	Product
CGS 8012.66	CareFRESH® Absorbent Sheets, 6.4 mm (1/4 in thick), 100 sheets
CGS 8013.66	CareFRESH® Absorbent Sheets, 3.2 mm (1/8 in thick), 200 sheets