

Wound Clips, Absorbable Sponges, and Bone Wax



Autoclip® Wound Clip System

This Autoclip® wound closure system provides a rapid, easy to use method to securely close wounds. When the wound is healed, simply remove the clip using the handy removal forcep. A convenient starter set is offered for first time buyers which includes an Autoclip® applicator, Autoclips®, and removal forcep. The stainless steel Autoclips® are available in two sizes, 9 mm and 18 mm. Instruments and wound clips are supplied non-sterile, but may be autoclaved.

- In use for over 40 years
- Provides maximum holding force with minimal trauma
- System includes 9 or 18 mm wound clips, applicators, and removal forcep
- Automatic, spring-fed applicators permit rapid, unassisted application

Catalog No.	Product
CGS 8858.67	Autoclip® Applicator for 9 mm Clips
CGS 8859.67	Autoclip® Applicator for 18 mm Clips
CGS 8860.67	Autoclips®, 9 mm, stainless steel, pkg. of 100
CGS 8861.67	Autoclips®, 18 mm, stainless steel, pkg. of 100
CGS 8862.67	Autoclip® Removal Forcep, for both 9 mm and 18 mm clips
CGS 8863.67	Autoclip® Starter Set, includes Applicator, 9 mm Autoclips, and Removal Forcep



- Highly absorbent
- Holds many times its weight in blood
- Completely absorbed within four to six weeks

Absorbable Gelatin Sponges

These small, sterile, surgical Sponges are prepared from specially treated, purified gelatin solution. Each Sponge is a hemostatic device that is capable of absorbing and holding within its meshes many times its weight in whole blood. These Sponges can be implanted during surgery and will be completely absorbed within four to six weeks.

These Gelatin Sponges measure 20 x 20 mm (3/4 x 3/4 in) and are 7 mm (1/4 in) thick. They are supplied sterile in a jar of 15 Sponges. They may be cut into smaller sizes.

Note: Gelatin Sponges should not be used in closure of skin incisions because they may interfere with healing of skin edges. These Sponges will expand so special care must be taken when placing in a cavity or in closed tissue spaces.

Catalog No.	Product
CGS 8864.67	Absorbable Gelatin Sponges



- To stop bleeding during surgery and to patch holes made in the skull

Bone Wax

Bone Wax is a sterile mixture of beeswax and isopropyl palmitate, a wax-softening agent, used to help control bleeding from bone surfaces. It can achieve local hemostasis of bone by acting as a mechanical (tamponade) barrier. It does not act biochemically and is minimally resorbable. Supplied in two sterile 2.5 g (0.09 oz) packets.

Catalog No.	Product
CGS 8865.67	Bone Wax