


## Disposable Electrodes



All of these Disposable Biopotential Electrodes exceed the AAMI (American Association for Medical Instrumentation) electrical requirements.

### Small Rewettable Electrodes

These 25 mm (1 in) diameter, disposable, silver/silver-chloride Electrodes are ideal for multiple use on small subjects or subjects with sensitive skin. Adhesion can be increased during application and decreased for removal by wetting the Electrode. The non-messy, bacteriostatic, conductive adhesive gel provides excellent electrical characteristics without irritating skin and maintains skin contact during patient movement. The gel entirely coats the subject side of the Electrode for secure adhesion.

These Electrodes provide a strong signal, stable baseline and superior electrical performance that exceeds all AAMI (American Association for Medical Instrumentation) electrical requirements. They are supplied as a package of 100 pouches each of which contains 3 electrodes on a vinyl carrier for a total of 300 electrodes. They attach to a gripper clip or standard snap connector leadwire.

- Gentle adhesive gel can be removed and reapplied multiple times for long term use
- For use with any ECG or bioelectric input lead cable with gripper clip or snap connectors

### Electrodes for Resting Subjects

The entire surface of these disposable Electrodes is coated with silver ink and an adhesive gel. This combination provides a cost effective Electrode with the performance of a standard silver/silver-chloride monitoring electrode. These Electrodes have rapid defibrillation recovery, provide stable baseline and have excellent electrical characteristics that exceed all AAMI (American Association for Medical Instrumentation) electrical requirements. The adhesive gel entirely coats the subject side of the Electrode providing secure attachment. This superior skin contact provides quality traces.

These 35 mm long x 22 mm wide (1.4 x 0.9 in) electrodes are supplied as a package of 10 vinyl carrying cards of 10 electrodes for a total of 100 electrodes.

- For use with any ECG or bioelectric input lead cable with banana plug connector (requires use of an alligator clip)

Catalog No.	Product
CGS 8712.71	Small Rewettable Electrodes

Catalog No.	Product
CGS 8713.71	Electrodes for Resting Subjects
CGS 8714.71	Alligator Clip, Connects to any 4 mm Banana Plug, pkg. of 10
CGS 8715.71	Preparatory Pads, individually packaged, 5.1 x 2.5 cm (2 x 1 in) Textured Cloth Pads soaked in 70% Isopropanol Solution with Pumice, prepares skin surface to accept Electrodes, pkg. of 1000