

Rocker II and Vacuum Blotter



Rocker II

The Rocker II performs gel staining and destaining, Northern, Southern and Western blots, mixes test tubes and rocks agglutination cards.

The Rocker II has an adjustable tilt angle 0 to 20 degrees. The speed is adjustable from 3 to 25 rpm. The load capacity is 2.7 kg (6 lbs). The ambient temperature range of the Rocker II is 4° to 65°C. A wall mount transformer supplies a low voltage power.

The Rocker II accepts a wide variety of flasks, plates and tubes. The tilt angle and speed can be adjusted without the need for tools. Another advantage is a thin power cord which easily passes through incubator door gaskets.

Comes supplied with a flat non-skid mat and an optional stacking tray and dimpled mat are available.

- Ideal for gel staining and destaining
- Northern, Southern, and Western blots
- Quiet and rugged
- Versatile
- Easy to use
- Excellent value
- Two-year warranty

Vacuum Blotter

The Vacuum Blotter performs a rapid transfer of DNA or RNA from an agarose gel on to suitable membrane.

It has a built-in vacuum pump and manometer ensure reproducible results without gel collapse. The pump is protected by a 3 liter buffer reservoir. It's simple to use clamping mechanism and easily manipulated plastic masks ensure success on the first try.

With an adjustable, precise control vacuum perfect results are ensured. The transfer efficiency is nearly 100% and the self contained unit is very compact yet still offers the largest transfer area. The vacuum Blotter is quiet during operation and easy to clean and maintain.

It comes supplied with two plastic masks and is also available in packages of five.

- Transfer nucleic acids in less than one hour
- Compact and easy to use
- Two-year warranty

Specifications

Dimension, H x W x D:	
Without Platform	12.7 x 19.1 x 17.8 cm (5 x 7.5 x 7 in)
With Platform	17.8 x 30.5 x 30.5 cm (7 x 12 x 12 in)
Power	115 VAC or 230 VAC input to 2 A, 12 VDC output adapters
Net Weight	3.9 kg (8.5 lbs)
Shipping Weight	5.9 kg (13 lbs)
Stacking Tray Clearance	6.3 cm (2.5 in)
Stacking Tray Weight	0.7 kg (1.5 lbs)

Catalog No.

Product

CGS 8375.39	Rocker II, 115 VAC, 60 Hz, with Rubber Non-Skid Mat
CGS 8376.39	Rocker II, 230 VAC, 50 Hz, CE, with Rubber Non-Skid Mat
CGS 8377.39	Stacking Tray with Rubber Non-Skid Mat
CGS 8378.39	Dimpled Mat

Specifications

Dimension, H x W x D	10 x 30 x 41 cm (4 x 12 x 16 in)
Power	115 VAC, 60 Hz, 15 W
Net Weight	4.5 kg (9.9 lbs)
Shipping Weight	6 kg (13.2 lbs)
Porous Transfer Area	28 x 28 cm (11 x 11 in)
Manometer	0 to 100 mbar

Catalog No.

Product

CGS 8379.39	Vacuum Blotter, 115 VAC, 60 Hz, 15 W
CGS 8380.39	Plastic Masks, pkg. of 5